


Acros Organics 钯催化剂和膦配体

钯催化剂和膦配体以其高催化活性和高选择性广泛应用于多种反应类型中，如不对称氢化，C-C偶联反应和C-杂偶联反应等。以下是Acros Organics 提供的常用的钯催化剂和膦配体。


| 英文名称 | 中文名称 | 货号 | 规格 |
|--|------------------------|-----------|--------|
| Tetrakis(triphenylphosphine)palladium(0), 99% | 四(三苯基膦)钯(0), 99% | 202380010 | 1 g |
| | | 202380050 | 5 g |
| | | 202380250 | 25 g |
| Palladium(II) acetate, 47.5% Pd | 醋酸钯(II), 47.5% Pd | 195180010 | 1 g |
| | | 195180020 | 2 g |
| | | 195180100 | 10 g |
| Palladium(II) chloride, 59% Pd | 氯化钯(II), 59% Pd | 195200010 | 1 g |
| | | 195200050 | 5 g |
| | | 195200250 | 25 g |
| Palladium(II) chloride, 99,999% | 氯化钯(II), 99,999% | 369670010 | 1 g |
| | | 369670050 | 5 g |
| Bis(triphenylphosphine)palladium(II) chloride, 15% Pd | 双(三苯基膦)氯化钯(II), 15% Pd | 197320010 | 1 g |
| | | 197320050 | 5 g |
| Bis(triphenylphosphine)palladium(II) chloride, 98% | 双(三苯基膦)氯化钯(II), 98% | 299252500 | 250 mg |
| | | 299250025 | 2.5 g |
| | | 299250050 | 5 g |
| | | 299250250 | 25 g |
| Tris(dibenzylideneacetone)dipalladium(0), 97% | 三(二亚苄基丙酮)二钯(0), 97% | 318775000 | 500 mg |
| | | 318770050 | 5 g |
| Tris(dibenzylideneacetone)dipalladium-chloroform adduct, 97% | 三(二亚苄基丙酮)二钯-氯仿加合物, 97% | 369340010 | 1 g |
| Bis(dibenzylideneacetone)palladium | 双(二亚苄基丙酮)钯 | 291970010 | 1 g |
| | | 291970050 | 5 g |
| Bis(triphenylphosphine)palladium(II) acetate, 99% | 双(三苯基膦)醋酸钯(II), 99% | 209270010 | 1 g |
| | | 209270050 | 5 g |
| Bis(tri-t-butylphosphine)palladium(0), 98% | 双(三叔丁基膦)钯(0), 98% | 363501000 | 100 mg |
| | | 363500010 | 1 g |
| Bis(tricyclohexylphosphine)palladium(0) | 双(三环己基膦)钯(0) | 363511000 | 100 mg |
| Allylpalladium chloride, dimer, 98% | 烯丙基氯化钯二聚物, 98% | 206835000 | 500 mg |
| | | 206830050 | 5 g |
| Bis(benzonitrile)palladium(II) chloride, 99+% | 双(氰苯基)氯化钯, 99+% | 207902500 | 250 mg |
| | | 207900010 | 1 g |
| | | 207900100 | 10 g |
| Tri-tert-butylphosphine, 95% | 三叔丁基膦, 95% | 360890010 | 1 g |
| | | 360890050 | 5 g |
| Tri-m-tolylphosphine, 98% | 三-间-甲苯基膦, 98% | 317280050 | 5 g |
| Tri-p-tolylphosphine, 95% | 三-对-甲苯基膦, 95% | 422330010 | 1 g |
| | | 422330050 | 5 g |
| Tri-o-tolyl-phosphine, 99% | 三-邻-甲苯基膦, 99% | 422320010 | 1 g |
| | | 422320050 | 5 g |
| | | 422320250 | 25 g |
| | | | |
| Tris(4-chlorophenyl)phosphine, 97% | 三(4-氯苯基)膦, 97% | 422210010 | 1 g |
| | | 422210050 | 5 g |
| Tris(4-methoxyphenyl)phosphine, 95% | 三(4-甲氧基苯基)膦, 95% | 422240010 | 1 g |
| | | 422240100 | 10 g |
| Trimesitylphosphine, 97% | 三异丙叉丙酮基膦, 97% | 321130010 | 1 g |
| | | 321130050 | 5 g |
| Triisopropylphosphine, 98% | 三异丙基膦, 98% | 317330010 | 1 g |
| | | 317330050 | 5 g |
| Tris(2,6-dimethoxyphenyl)phosphine, 97+% | 三(2,6-二甲氧基苯基)膦, 97+% | 364960050 | 5 g |
| | | 364960250 | 25 g |
| Tricyclohexylphosphine, 97% | 三环己基膦, 97% | 421610010 | 1 g |
| | | 421610050 | 5 g |
| | | 421610250 | 25 g |
| Tri-(2-furyl)phosphine, 98% | 三(2-呋喃基)膦, 98% | 294802500 | 250 mg |
| | | 294800010 | 1 g |
| 1,1'-Bis-(diphenylphosphino)ferrocene | 1,1'-双(二苯基膦)二茂铁 | 348010050 | 5 g |
| | | 348010250 | 25 g |

Acros Organics 钯催化剂和磷配体


| 英文名称 | 中文名称 | 货号 | 规格 |
|---|------------------------------------|-----------|--------|
| 1,2-Bis-(diphenylphosphino)ethane, 98+% | 1,2-双-(二苯基膦基)乙烷, 98+% | 147910010 | 1 g |
| | | 147910100 | 10 g |
| | | 147910500 | 50 g |
| 1,3-Bis(diphenylphosphino)propane, 97% | 1,3-双(二苯基膦基)丙烷, 97% | 310050010 | 1 g |
| | | 310050050 | 5 g |
| | | 310050250 | 25 g |
| 1,4-Bis-(diphenylphosphino)butane, 98% | 1,4-双-(二苯基膦基)丁烷, 98% | 296460050 | 5 g |
| | | 296460250 | 25 g |
| 1,5-Bis(diphenylphosphino)pentane, 97% | 1,5-双(二苯基膦)戊烷, 97% | 320850010 | 1 g |
| | | 320850050 | 5 g |
| 2-(Di-tert-butylphosphino)biphenyl, 99% | 2-(二叔丁基膦基)联苯, 99% | 356215000 | 500 mg |
| | | 356210050 | 5 g |
| | | 356210250 | 25 g |
| 2-(Dicyclohexylphosphino)biphenyl, 98% | 2-(二环己基膦基)联苯, 98% | 356225000 | 500 mg |
| | | 356220020 | 2 g |
| | | 356220050 | 5 g |
| 2-(Dicyclohexylphosphino)-2'-(N,N-dimethylamino)biphenyl, 97% | 2-(二环己基膦基)-2'-(N,N-二甲胺基)联苯, 97% | 356235000 | 500 mg |
| | | 356230010 | 1 g |
| | | 356230050 | 5 g |
| 2-(Dicyclohexylphosphino)-2',4',6'-trisopropylbiphenyl, 98+% | 2-(二环己基膦基)-2',4',6'-三异丙基联苯, 98+% | 387140010 | 1 g |
| | | 387140050 | 5 g |
| | | 387140250 | 25 g |
| (R)-(+)-2,2'-Bis(diphenylphosphino)-1,1'-binaphthyl, 99+% | (R)-(+)-2,2'-双(二苯基膦)-1,1'-联萘, 99+% | 265530250 | 25 mg |
| | | 265532500 | 250 mg |
| | | 265530010 | 1 g |
| | | 265530050 | 5 g |
| (S)-(-)-2,2'-Bis(diphenylphosphino)-1,1'-binaphthyl, 99+% | (S)-(-)-2,2'-双(二苯基膦)-1,1'-联萘, 99+% | 265542500 | 250 mg |
| | | 265540010 | 1 g |
| | | 265540050 | 5 g |
| (±)-2,2'-Bis(diphenylphosphino)-1,1'-binaphthyl, 98% | (±)-2,2'-双(二苯基膦)-1,1'-联萘, 98% | 368640010 | 1 g |
| | | 368640050 | 5 g |
| | | 368640250 | 25 g |
| Tri-tert-butylphosphonium tetrafluoroborate, 97+% | 三叔丁基膦四氟硼酸, 97+% | 366940010 | 1 g |
| | | 366940050 | 5 g |




20238




20790



36089




34801




14791




31005




29480




42232




36864



35621



42161



42221



Fisher Scientific

Part of Thermo Fisher Scientific

For customer service, call 400-881-5117.
To send Email: sales.china@thermofisher.com
To order online: www.fishersci.com.cn

©2010 Thermo Fisher Scientific Inc. All rights reserved.